


alpress[®]
technology for health



IDEA
AERO DUO

Respiratory Ventilators



alpress®
technology for health



IDEEA
AERO DUO


Respiratory Ventilators

EON
GROUP



Ventilator ICU Ventilation and Respiratory Monitoring

Fully-featured, high-performance ventilation with modern Invasive and Non-invasive ventilation modes



IDEA
AERO DUO
Respiratory Ventilators

EON
GROUP

100%

Made In Türkiye

IDEA AERO DUO



IDEA Aero is produced in TURKEY by Alpress with the support of biomedical engineers and ICU experienced medical doctors
Aero is the result of 17 years high technology machines manufacturing experience.

AERO DUO HIGH LEVEL INTENSIVE CARE MECHANICAL VENTILATOR OF TURKEY



DESIGNED AND PRODUCED FOR INTENSIVE CARE AND
REANIMATION UNITS.

- Modern and ergonomic lines.
- User-friendly interface.
- Perfect performance.
- Traditional and innovative operating modes.
- Compatible with pediatric and adult patients.
- Low cost and maintenance.

IDEA
AERO DUO
Respiratory Ventilators

EON
GROUP



Ventilator ICU Ventilation and Respiratory Monitoring

Monitoring Functions

- Real time FiO2 monitoring
- Volume, Flow and Pressure waveform graphics
- Leakage monitoring and compensation
- Up to 8 hours battery life
- Wifi support for online monitoring
- Peep Level monitoring
- Plateau reading

Key Features

- Non-invasive ventilation in all modes
- 2-30 l/min %100 O2 Therapy without time limit
- Full record of all ventilation data
- Easy to use features with low-learning curve
- APRV Ventilation mode

Advanced Safety

- Detachable exhalation valve for avoid cross-contamination
- 4 flow sensor for advanced accuracy
- Spontaneous breath triggering with pressure and flow sensors
- 30 - 2500ml tidal volume for every patient needs
- Intelligent alarms for maximum patient safety
- Auto calibration and diagnostic for stability
- Up to 8 hours battery life




Respiratory Ventilators



Ventilator ICU Ventilation and Respiratory Monitoring

Pediatric
and Adult

One Machine Compatible with All Patients



1. Alarm limits
2. Active ventilation mode
3. Realtime respiratory values
4. Message and status bar
5. Battery level
6. Date & Time
7. Alarm silent button
8. Sub-menu buttons
9. Fast ventilation values see / modify area
10. Pressure waveform
11. Flow waveform
12. Volume waveform

Advanced
Adaptive Control

Aero works both in traditional and
innovative modes.



Pressure Controlled
Ventilation

- PC-CMV
- PC-SIMV
- PC-AC
- PC-APRV
- PV-PSV

Volume Controlled
Ventilation

- VC-CMV
- VC-SIMV
- VC-AC

Support of Spontaneous
Breathing

- SPN-CPAP
- SPN-BILEVEL CPAP

Advanced
Adaptive Control

Aero uses advanced adaptive control algorithms.
It responds in milliseconds.





Ventilator ICU Ventilation and Respiratory Monitoring




International Standard Flanges

Easy connect air and oxygen flanges doesn't require any tools and fits any wall plug connections.



Removable and Autoclavable Patient Circuit Connections

For avoiding cross-contamination of patients Aero duo circuit connections removable and disinfectable.



High Resolution Touch Screen

10.1 inch high resolution lcd touch screen of IDEA duo presents clear and highly visible graphics and respiratory values from any view angle.

Also monitor can be rotate up to users position.



Rotary Button for Easy Adjustment

High sensitivity and robust rotary button provides accurate adjustment of all values of respiration modes. Also push button feature makes confirmation easy.



Ventilator ICU Ventilation and Respiratory Monitoring



Pediatric	Adult
Respiratory rate: 1-120 Per Minute	Respiratory rate: 1-80 Per Minute
Tidal Volume: 0.02-06 Ltr	Tidal Volume: 0.02-06 Ltr
Flow: 1-80 ltr per minute	Flow: 1-120 ltr per minute
Trigger Sensitivity: 0.1-20 ltr per minute 0.1-20 mbar	Trigger Sensitivity: 0.1-20 ltr per minute 0.1-20 mbar
Leak Compensation	Leak Compensation
Invasive and Non-Invasive Ventilation	Invasive and Non-Invasive Ventilation



Ventilator ICU Ventilation and Respiratory Monitoring

Technical Features of Aero

Types of Patients:

Padiatric and Adult



Modes of Ventilation

- VC-CMV
- VC-SIMV
- VC-AC
- PC-CMV
- PC-SIMV
- PC-AC
- PC-APRV
- PV-PSV
- SPN-CPAP
- SPN-BILEVEL CPAP

Detailed Features

- Inspiratory Time (Thigh) • 0.1 to 30 s
- Expiratory Time (Tlow) • 0.05 to 30 s
- Inspiratory Pressure (Phigh) • 1 to 95 mbar (or hPa or cmH₂O)
- Expiratory Pressure (Plow) • 0 to 50 mbar (or hPa or cmH₂O)
- Airway Presure Measurement
 - Auto Calibration and Auto Diagnostic
 - Realtime Flow Calibration
 - Auto Leakage Compensation
 - 2-30l/min O₂ Therapy
 - Both invasive and non invasive modes
- Enhancements
 - Auto Calibration and Auto Diagnostic
 - Realtime Flow Calibration
 - Auto Leakage Compensation
 - 2-30l/min O₂ Therapy
 - Both invasive and non invasive modes

Working Features

- Therapy Types
 - Invasive Ventilation (Tube)
 - Non-invasive ventilation (NIV)
 - O₂-Therap0y
- Ventilation Frequency (RR)
 - Adult 0.5 to 80 / min
 - Pediatric patients, 0.5 to 80 / min
- Inspiration Time (Ti)
 - Adults 0.2 to 10 s
 - Pediatric Patients, 0.2 to 10 s
- Tidal Volume (VT)
 - Adults 0.1 to 2.5L
 - Pediatric Patients 0.03 to 0.3L
- Inspiratory Flow (Flow)
 - Adults 2 to 120L / min
 - Pediatric Patients, 2 to 30L / min
- Inspiratory Pressure (Pinsp) • 1 to 95 mbar (or hPa or cmH₂O)
- Inspiratory Pressure Limit (Pmax) • 2 to 100 mbar (or hPa or cmH₂O)
- PEEP • 0 to 50 mbar (or hPa or cmH₂O)
- Pressure Assist (Psupp) • 0 to 95 mbar (or hPa or cmH₂O)
- Rise time for pressure assist • Adults, Pediatric patients 0 to 2s
- O₂ Concentration (FIO₂) • 21 to 100 Vol. %
- Triger sensitivity (Flow Trigger) • 0.2 to 15L / min



Ventilator ICU Ventilation and Respiratory Monitoring

Flow Measurement

- | | |
|-----------------------------------|--|
| Minute Volume Measurement | <ul style="list-style-type: none"> • Expiratory Minute Volume (MVe) • Inspiratory Minute Volume (MVi) • Total Minute Volume (MV) |
| Tidal Volume Measurement | <ul style="list-style-type: none"> • Tidal Volume (VT) • Inspiratory Tidal Volume • Expiratory Tidal Volume |
| Respiratory Rate Measurement | <ul style="list-style-type: none"> • Breathing Frequency (RR) • Mandatory Respiratory Rate (RRmand) • Spontaneous Breathing Frequency (RRspont) |
| O2 Measurement (Inspiratory Side) | <ul style="list-style-type: none"> • Inspiratory O2 Concentration (Fio2) • Range 18 to 100 Vol% |

- Curve Displays
- Airway Pressure Paw (l) -30 to 100 mbar (or hPa or cmH2O)
 - Flow (l) -100 to 100 L / min
 - Volume V (l)

Mains Power Supply

Mains power connection • 100 V to 240 V, 50/60 Hz

Current Consumption

At 230 V • Max. 0.5 A
At 100 V • Max. 0.9 A

Alarms / Monitoring

- | | |
|---|-------------------|
| Expiratory Minute Volume (MVe) | • High / Low |
| Airway Pressure (Paw) | • High / Low |
| Inspiratory O2 Concentration (Fio2) | • High / Low |
| Volume Monitoring (VT) | • High / Low |
| Apnea Alarm Time (Tapn) | • 5 to 60 seconds |
| Disconnect Alarm Delay Time (Tdisconnect) | • 0 to 60 seconds |

Performance Data

- | | |
|-------------------------------|--|
| Control Principle | • Time-cycled, Volume-constant, Pressure-controlled |
| Inspiratory Flow | • Max, 180 L / min |
| Base flow, adults | • 2 L / min |
| Base flow, pediatric patients | • 3 L / min |
| Safety Valve | • Opens if medical compressed air supply fails (supply gas flow is not sufficient to provide the inspiratory flow required), enables spontaneous breathing with ambient air. |

Power Consumption

- | | |
|--|---|
| Maximum | • 100 W |
| During Ventilation, Without Charging the Battery | • Approx. 80 W |
| Digital Machine Output | • Digital output and input via an RS232 C interface |

Gas Supply

- | | |
|--------------------|--------------------------------------|
| O2 Gauge Pressure | • 2.0 to 6.0 bar (or 200 to 600 kPa) |
| Air Gauge Pressure | • 2.0 to 6.0 bar (or 200 to 600 kPa) |

Size and Weight


- | | |
|-----------------------|----------------------------|
| Aero Duo | • Approx. 16 kg (35.2 lbs) |
| Aero Duo With Trolley | • Approx. 34 kg (74.9 lbs) |

Physical Specifications

- | | |
|---------------------------------|--|
| Ventilation Unit | • 460 mm x 210 mm x 310 mm
(18.1 in x 8.2 in x 12.2 in) |
| Ventilation Unit on the Trolley | • 500 mm x 1210 mm x 400 mm
(19.6 in x 47.6 in x 15.7 in) |

Aero Duo Screen

Diagonal Screen Size • 10.1" TFT Color Touch Screen



Diagonal Screen Size • 10.1" TFT Color Touch Screen



Aero Duo Screen

Diagonal Screen Size • 10.1" TFT Color Touch Screen

IDEA
AERO DUO
Respiratory Ventilators



Ventilator ICU Ventilation and Respiratory Monitoring



Comparison Of Mode



IDEA AERO	Puritan Bennett	Dräger	GE	Hamilton	Maquet	Mindray	Aeonmed
VCV	A/C:PC	PC-AC	PCV	PCV	PC	P-AC	VCV-AC
PSV	SIMV: PC	PC-SIMV	SIMV-PC	PSIMV+	SIMV-PC+PS	P-SIMV	PSV-AC
AC	PS	SPN-CPAP/PS	CPAP/PSV	Spont	PS	PSV	SIMV/VCV+PSV
SIMV / VC	BiLevel	PC-BIPAP	BiLevel	DuoPAP	Bi Vent	DuoLevel	SIMV/PCV+PSV
SIMV / PS	APRV	PC-APRV	APRV	APRV	Bivent-APRV	APRV	SIMV/PRVC+PSV
SIMV / PSV	A/C: VC	VC-AC	VCV	(S)CMV	VC	V-AC	SPONT/CPAP+PSV
CPAP	VC+	Autoflow	PCV-VG	APV/SIMV+	PRVC	PRVC	PSV / VG
Bilevel CPAP	A/C: VC	VC-CMV	VCV	CMV	VC	V-AC	SPON / CPAP + PSV
Hi Level O2 Therapy	SIMV: VC	VC-SIMV	SIMV-VC	SIMV	SIMV-VC+PS	V-SIMV	
PSV / VG	VC+	VC-SIMV+ Autoflow	SIMV-PCVG	APV/SIMV+	SIMV-PRVC+PS	PRVC	
-	PS	SPN-CPAP/PS	CPAP	Spont.	PS/CPAP	-	
-	VS	SPN-CPAP/VS	-	-	VS	-	
-	PAV+	SPN-PPS	-	-	-	-	
-	-	-	-	ASV	-	-	
-	-	-	-	-	NAVA	-	



Ventilator ICU Ventilation and Respiratory Monitoring

FRONT SIDE CONNECTING POINTS




Figure 2: Connections at the front side of the device

- 1 **Connection of tube circuit - inspiration**
The single line patient circuit or the inspiration section of a double line patient circuit is connected here.
- 2 **Connection of pressure measuring tube**
The single line patient circuit or the inspiration section of a double line patient circuit is connected here.
- 3 **Connection of tube circuit - expiration**
The expiration section of a double line patient circuit will be connected here.
- 4 **Expiration flow sensor input socket**
Expiration sensor cable will connect here

REAR SIDE CONNECTING POINTS




Figure 3: Connections at the rear side of the device

- 5 **Oxygen sensor connection**
- 6 **Electrical Connection 110-220V 0.8Amax**
- 7 **Medical Oxygen input Max 6 Bar**
- 8 **Medical Air input Max 6 Bar**

8 Medical Air input Max 6 Bar


7 Medical Oxygen input Max 6 Bar



Ventilator ICU Ventilation and Respiratory Monitoring

OPTIONAL BLOWER / TRIBUNE DESCRIPTION

IDEA AERO BLOWER DRIVE



- 1- Pressurized air outlet
- 2- Air inlet filter group
- 3- Blower electronics and sensor housing
- 4- Power input and hall sensor output connection

CONTROL ELEMENTS



Figure 4: Control elements

- 9 3 color Led indicator
- 10 Rotary encoder and push button

10 Rotary encoder and push button

9 3 color Led indicator



Ventilator ICU Ventilation and Respiratory Monitoring

EXPIRATION VALVE SET



Expiration

- 11 Set Detachable Flow sensor
- 12 Disposable-Reusable expiration valve

MOVABLE AND REMOVABLE HOUSING PARTS




Figure 6: Device Handle

- 14 Handle (pull-out)
The handle may be pulled out for device transport.




Figure 7 :Control Monitor


- 15-Control Monitor,
Monitor is rotatable



Ventilator ICU Ventilation and Respiratory Monitoring

CONNECTING A DOUBLE LINE PATIENT CIRCUIT

Connect the tube circuit to the device according



A Patient side connection B Pressure measuring tube C Expiration tube D Inspiration tube E+H Bacterial filters


Figure 16: Connecting a double line patient circuit with pressure measuring tube

FiO₂ SENSOR CONNECTION

Connect the FiO₂ sensor to the device as shown in Figure 23. If ventilation is in progress, it will also show in the parameter screen.

NOTICE

In the case of a malfunction (abnormal operation) dashes are displayed instead of the measured value.



A Oxygen Sensor
B Oxygen Sensor Cable
C Oxygen Sensor Socket
D Oxygen Cable Socket

D Oxygen cable socket
C Oxygen sensor socket
B Oxygen sensor cable
A Oxygen sensor




Ventilator ICU Ventilation and Respiratory Monitoring

CONNECTING OXYGEN AND AIR SOURCE

ATTENTION

Only the oxygen connection adapter supplied may be used to connect oxygen. Otherwise, there is a risk that the back-stop in the connection is damaged.

Connect the oxygen source to the device as shown in Figure 27.



A Tube From oxygen source B Tube from Medical Air source, C Medical Air Connection D Oxygen connection

SWITCHING THE DEVICE ON

NOTICE

The tube circuit may be connected when the device is started up, but it may not yet be connected to the patient yet.

If you are using oxygen therapy during ventilation, please note the section "Using oxygen"

To switch on the device:

1. Press the main power switch on the back of the device (position "ON").



Figure 31: Switching on the device

2. The home screen will be displayed.

SWITCHING THE DEVICE OFF

1. Stop the ventilation.
2. Switch off the power with the main power switch on the rear panel (position "OFF").




Figure 32: Switching off the device


©2021 ALPRESS All rights reserved.

The content published in this user's manual is the sole property of the ALPRESS . All forms of reproduction, editing, distribution and any kind of exploitation, even in part, require the prior written consent of the ALPRESS .

The ALPRESS reserves the right to amend or replace this user's manual without prior notice.

Please ensure that you are always working with the most current version of this user's manual. Should you have any questions, please contact the ventilation device provider, or check our information at www.eon-group.com

The respiratory device may only be operated and maintained by trained personnel.



About Us

EON GROUP are global distributor of IDEA AREO DUO on behalf of ALPRESS company. We are global distributor, for more detail please visit our group website | www.eon-group.com

alpress

TASARIM KALIP DANIŞMANLIK LTD. ŞTİ

Seyit Nizam Mah. Demirciler Sitesi 9.Yol No:26
Zeytinburnu/ İstanbul/ TURKEY
Tel: +90 212 416 65 05 - 416 65 08
Fax: +90 212 416 65 54
info@alpress.com.tr | www.alpress.com.tr

To whom it may concern,

We, **ALPRESS DESIGN MOULD & CONSULTING CO.** hereby declare that **EON GROUP (THAILAND) CO., LTD.** is our Exclusives Authorization Dealer to the following countries.

Thailand, India, and all over Asia

This Authorization letter can be verified by our head office in ISTANBUL.
Sincerely,

18.08.2021

Alpress Kalıpcılık Dan. İth. İhr. Ltd. Şti.


KALIPÇILIK DANIŞ. İTH. İHR. SAN. VE TİC. LTD. ŞTİ.
Demirciler Sit. 9. Yol No: 26 Zeytinburnu / İST.
Tic. Sic. No: 275801/08 / 08 / Fax: 0212 416 65 54
Davutpaşa / İD: 258 0127 032


Our License

EON GROUP are global distributor of
IDEA AREO DUO on behalf of ALPRESS
company. We are global distributor, for
more detail please visit our group
website | www.eon-group.com




Unique

High Quality Machine
made by ALPRESS




First to Market

First one manufacture
100% Made in Turkey



Tested

Certify under
ISO13485, CE, FDA



Authentic

SARS-COV-2, invasive
ventilation is required.



Product


IDEA AERO DUO Ventilator is
ICU Ventilation and Respiratory
Monitoring machine made by
ALPRESS Company.

100% made in Turkey.

IDEA AERO DUO

ICU VENTILATOR AND RESPIRATORY

Respiratory Ventilators



FAST PRODUCED VENTILATOR SYSTEM

As alpress,

We have produced a minimal (clinically accepted) ventilator production in hospitals related to the current pandemic caused by the COVID-19 virus.

Aero Duo for adults and children from a tidal volume of 50 ml and higher and in professional healthcare facilities when an FDA-cleared clinical ventilator is not available during the COVID-19 pandemic. It can be used for intensive care ventilation. The Aero Duo is not suitable for use in vehicles, airplanes and helicopters.

Alpress reserves the right to change any products, technical specifications, price and stock information, without any prior notice.

Alpress, istediği zaman ve önceden herhangi bir bildirimle gerek olmaksızın modelleri, donanımı, teknik özellikleri, fiyat ve stok bilgilerini değiştirme hakkını saklı tutar.



IDEA AERO
Respiratory Ventilators

Technology For Health




IDEA Aero Respiratory Ventilators

Alpress produced the IDEA Aero Respiratory Ventilators by proving itself once again with its technological infrastructure and knowledge





IDEA AERO
Respiratory Ventilators



Ventilation Modes

- VCV
- PCV
- AC-VC/PC
- SIMV
- SIMV/PC
- CPAP



Extensive monitoring



- 3 waveforms
- Alarm History
- 30 days patient trend data
- Up to 5000 events log



Technology For Health



IDEA AERO
Respiratory Ventilators



IDEA Aero comes with adjustable screen that allows tilting for best viewing angle



- Setting range indicator
- Setting summary



- Alarm range indicator
- Active alarm with color coded
- Auto alarm limits



Easy of use

Thanks to the innovative UI design of IDEA Aero, change of ventilation modes requires only 2 simple steps. Each function is in logical order so that you wouldn't be lost in complex user manual.



Technology For Health

Ease of maintenance

With our newly designed expiratory and inspiration valves, replacing and cleaning are no longer a hassle. Also, it is extremely durable and can last longer than most of the competitors offerings right now out there in the market.

Valve system:

- Detachable design
- Require no tools during disassemble
- Support autoclave disinfection
- Extremely durable






IDEA AERO
Respiratory Ventilators

Technology For Health


We know that uninterrupted patient care is important. With our IDEA Aero mechanical ICU ventilator, your patients are always under support and protection no matter how stressful the environments are. Also, with our powerful patient monitoring solution, patient data can be retained and transfer seamlessly.




- Easy detachable trolley
- Light weight only ~16 Kg
- Up to 4 hrs internal battery support
- FIO2 monitoring

- Modern ventilation modes
- Waveforms - Trends
- Advanced Safety
- Long term Logging

- Multi touch user interface
- Silent operation
- Better synchronous breaths
- Detachable exhalation valve



Aero Duo



IDEA Aero Uses Simple but Powerfull Interface to Ensure Stability

- IDEA Aero uses 7 most common control modes.
- Volume and pressure controlled modes can be selected on touch screen
- O2 - Air mix can be adjusted from touch screen and O2 concentration can be adjusted 21%-100%
- IDEA Aero works compitable with all adult patients.
- 50-3000 ml Vtial and up to 80 breaths in a minute
- Flow and Pressure graphs to inspect ventilation process



IDEA AERO
Respiratory Ventilators



Aero Duo

Monitoring Functions

- Real time FiO₂ monitoring
- Volume, Flow and Pressure waveform graphics
- Leakage monitoring and compensation
- Up to 8 hours battery life (optional)
- Peep Level monitoring
- Plateau reading

Advanced Safety

- Detachable exhalation valve for avoid cross-contamination
- 4 flow sensor for advanced accuracy
- Spontaneous breath triggering with pressure and flow sensors
- 50 - 3000ml tidal volume for every patient needs
- Intelligent alarms for maximum patient safety
- Auto calibration and diagnostic for stability

Key Features

- Non-invasive ventilation in all modes
- Full record of all ventilation data
- Easy to use features with low-learning curve
- Multi language





Ventilation Setting

Ventilation Mode		
Volume Controlled Ventilation	Pressure Controlled Ventilation	Support of Spontaneous Breathing
<ul style="list-style-type: none">• VCV• AC/VC• SIMV-VC	<ul style="list-style-type: none">• PCV• AC/PC• SIMV/PC	<ul style="list-style-type: none">• SPN-CPAP
Enhancements	Therapy Types	
<ul style="list-style-type: none">• Auto Calibration and Auto Diagnostic• Realtime Flow Calibration• Auto Leakage Compensation• Both invasive and non invasive modes	<ul style="list-style-type: none">• Invasive Ventilation (Tube)• Non-invasive ventilation (NIV)	
Ventilation Frequency (RR)	Inspiration Time (Ti)	
<ul style="list-style-type: none">• Adult 0.5 to 80 / min• Pediatric patients, 0.5 to 80 / min	<ul style="list-style-type: none">• Adults 0.2 to 10 s• Pediatric Patients, 0.2 to 10 s	
Tidal Volume (VT)	Inspiratory Flow (Flow)	
<ul style="list-style-type: none">• Adults 0.1 to 2.5L• Pediatric Patients 0.03 to 0.3L	<ul style="list-style-type: none">• Adults 2 to 120L / min• Pediatric Patients, 2 to 30L / min	
Inspiratory Pressure (Pinsp)	Inspiratory Pressure Limit (Pmax)	
<ul style="list-style-type: none">• 1 to 95 mbar (or hPa or cmH2O)	<ul style="list-style-type: none">• 2 to 100 mbar (or hPa or cmH2O)	
PEEP	Pressure Assist (Psupp)	
<ul style="list-style-type: none">• 0 to 50 mbar (or hPa or cmH2O)	<ul style="list-style-type: none">• 0 to 95 mbar (or hPa or cmH2O)	
Rise time for pressure assist	O2 Concentration (FiO2)	
<ul style="list-style-type: none">• Adults, Pediatric patients 0 to 2s	<ul style="list-style-type: none">• 21 to 100 Vol. %	
Triger sensitivity (Flow Trigger)		
<ul style="list-style-type: none">• 0.2 to 15L / min		

Inspiratory Time (Thigh)	Expiratory Time (Tlow)
<ul style="list-style-type: none">• 0.1 to 30 s	<ul style="list-style-type: none">• 0.05 to 30 s
Inspiratory Pressure (Phigh)	Expiratory Pressure (Plow)
<ul style="list-style-type: none">• 1 to 95 mbar (or hPa or cmh2O)	<ul style="list-style-type: none">• 0 to 50 mbar (or hPa or cmh2O)

Displayed Measured Values

Airway Presure Measurement
<ul style="list-style-type: none">• Plateau Presure (Pplat)• Positive end-expiratory Pressure (PEEP)• Positive end-expiratory Pressure (PIP)• Mean airway pressure (Pmean)• Range -30 to 100 mbar (or hPa or cmH2O)

Flow Measurement

Minute Volume Measurement	Tidal Volume Measurement
<ul style="list-style-type: none">• Expiratory Minute Volume (MVe)• Inspiratory Minute Volume (MVi)• Total Minute Volume (MV)	<ul style="list-style-type: none">• Tidal Volume (VT)• Inspiratory Tidal Volume• Expiratory Tidal Volume
Respiratory Rate Measurement	O2 Measurement (Inspiratory Side)
<ul style="list-style-type: none">• Breathing Frequency (RR)• Mandatory Respiratory Rate (RRmand)• Spontaneous Breathing Frequency (RRspon)	<ul style="list-style-type: none">• Inspiratory O2 Concentration (Fio2)• Range 12 To 100 Vol%
Curve Displays	
<ul style="list-style-type: none">• Airway Pressure Paw (t) -30 to 100 mbar (or hPa or cmH2O)• Flow (t) -100 to 100 L / min• Volume V (t)	




IDEA AERO
Respiratory Ventilators

Alarms / Monitoring

Expiratory Minute Volume (MVe) <ul style="list-style-type: none">• High / Low	Airway Pressure (Paw) <ul style="list-style-type: none">• High / Low
Inspiratory O2 Concentration (Fio2) <ul style="list-style-type: none">• High / Low	Volume Monitoring (VT) <ul style="list-style-type: none">• High / Low
Apnea Alarm Time (Tapn) <ul style="list-style-type: none">• 5 to 60 seconds	Disconnect Alarm Delay Time (Tdisconnect) <ul style="list-style-type: none">• 0 to 60 seconds

Performance Data

Control Principle <ul style="list-style-type: none">• Time-cycled, Volume-constant, Pressure-controlled	Inspiratory Flow <ul style="list-style-type: none">• Max. 120 L / min
Base flow, adults <ul style="list-style-type: none">• 2 L / min	Base flow, pediatric patients <ul style="list-style-type: none">• 3 L / min
Safety Valve <ul style="list-style-type: none">• Opens if medical compressed air supply fails (supply gas flow is not sufficient to provide the inspiratory flow required), enables spontaneous breathing with ambient air.	



Aero Duo

Operating Data

Mains Power Supply

Mains power connection

- 100 V to 240 V, 50/60 Hz

Current Consumption

At 230 V	At 100 V
<ul style="list-style-type: none">• Max. 0.5 A	<ul style="list-style-type: none">• Max. 0.9 A

Power Consumption

Maximum	During Ventilation, Without Charging the Battery
<ul style="list-style-type: none">• 100 W	<ul style="list-style-type: none">• Approx. 80 W
Digital Machine Output <ul style="list-style-type: none">• Digital output and input via an RS232 C interface	

Gas Supply

O2 Gauge Pressure	Air Gauge Pressure
<ul style="list-style-type: none">• 2.0 to 6.0 bar (or 200 to 600 kPa)	<ul style="list-style-type: none">• 2.0 to 6.0 bar (or 200 to 600 kPa)

Physical Specifications

Ventilation Unit	Dimensions (W x H x D) <ul style="list-style-type: none">• 460 mm x 210 mm x 310 mm (18.1 in x 8.2 in x 12.2 in)
Ventilation Unit on the Trolley	<ul style="list-style-type: none">• 500 mm x 1210 mm x 400 mm (19.6 in x 47.6 in x 15.7 in)

Weight

Aero Duo	<ul style="list-style-type: none">• Approx. 16 kg (35.2 lbs)
Aero Duo With Trolley	<ul style="list-style-type: none">• Approx. 34 kg (74.9 lbs)

Aero Duo

Diagonal Screen Size	<ul style="list-style-type: none">• 10.1" TFT Color Touch Screen
----------------------	--



IDEA AERO DUO²

OXYGEN CONCENTRATOR

Oxygen Concentrator



SECOND PRODUCED OXYGEN CONCENTRATOR

As alpress,

We have produced a minimal (clinically accepted) ventilator production in hospitals related to the current pandemic caused by the COVID-19 virus.

Aero Duo for adults and children from a tidal volume of 50 ml and higher and in professional healthcare facilities when an FDA-cleared clinical ventilator is not available during the COVID-19 pandemic. It can be used for intensive care ventilation. The Aero Duo is not suitable for use in vehicles, airplanes and helicopters.

Alpress reserves the right to change any products, technical specifications, price and stock information, without any prior notice.

Alpress, istediği zaman ve önceden herhangi bir bildirme gerek olmaksızın modelleri, donanımı, teknik özellikleri, fiyat ve stok bilgilerini değiştirme hakkını saklı tutar



IDEA AERO²
Oxygen Concentrator

Technology For Health



IDEA Aero² Oxygen Concentrator

Alpress produced the IDEA Aero² Oxygen Concentrator by provroving itself once again with its technological infrastructure and knowledge






IDEA AERO²
Oxygen Concentrator



Oxygen Concentrator


- 10 LPM Type
- 5 LPM Type
- 5.5 inch touchscreen
- 24 hour continuous
- Medical grade
- Humidifier



Technology For Health



IDEA AERO²
Oxygen Concentrator



Easy of use

Thanks to the innovative UI design of IDEA Aero Duo 2, Change of Oxygen Concentrator control only 2 simple steps, By 5.5 inch touchscreen control is in logical order so that you wouldn't be lost in complex user manual.



Technology For Health

Ease of maintenance

With our newly designed expiratory and inspiration valves , replacing and cleaning are no longer a hassle. Also, it is extremely durable and can last longer than most of the competitors offerings right now out there in the market.

Valve system:

- Detachable design
- Require no tools during disassemble



IDEA AERO²
Oxygen Concentrator



Aero Duo

Product Type

- 10 LPM
- 5 LPM

Liter flow

- 0.5 to 10 liters per minute high purity medical grade O2 supply
- 0.5 to 5 liters per minute high purity medical grade O2 supply

Product Specifications

- 24 hour continuous work capacity
- Medical grade oil free silent compressor
- Built-in Humidifier
- External compressor connection option
- Shock absorber lockable wheels
- 5.5 inch touchscreen control

Flow Measurement


- 42 dBA typical
 - 40 dBA typical
- OPI (oxygen percentage indicator) alarm levels
- Low oxygen: 82%, very low oxygen: 70%
- Operating temperature
- 55°F to 90°F (12°C to 32°C)
- Storage/transport humidity
- -30°F to 160°F (-34°C to 71°C) up to 95% relative humidity
- Operating humidity
- Up to 95% relative humidity
- Outlet pressure
- 5.5 PSI
- Operating altitude
- 0 to 7,500' (0 to 2,286 m)

Power Consumption

- 220 VAC +/- 10% 50hz
- Lightweight (21Kg - 46lbs)
- Compact (45cmX20cmX45cm)
- 100% Made in Turkey
- 3 year guarantee
- Lightweight (17Kg - 46lbs)

Products Certification


Idea Aero Duo made under certification follow:



Products Certification

Idea Aero Duo made under certification follow:

Test Report




NETES MÜHENDİSLİK VE DİŞ TİC. A.Ş.
KALİBRASYON LABORATUVARI
Küçük Çamlıca Mah. Oymak Çıkmazı Sok.
No:3 Üsküdar / İstanbul

KALİBRASYON SERTİFİKASI
Calibration Certificate

B2000350
NETES KALİBRASYON LABORATUVARI
05-20


Cihaz Adı Device Name	: VENTİLATÖR CİHAZI
Markası (Üreticisi) Trademark (Manufacturer)	:
Tipi / Modeli Type / Model	: DV0104
Seri Numarası Serial Number	: ADMY2000250
Demirbaş Numarası Device ID	:
İstek No Order No	: 20-01189
Cihaz Sahibi (Adı / Adresi) Customer (Name / Address)	:
Kalibrasyon Tarihi Date of Calibration	: 27.5.2020
Sertifika Sayısı Number of pages of the Certificate	: 3

Bu kalibrasyon sertifikası, Uluslararası Birimler Sisteminde (SI) tanımlanmış birimleri realize eden ulusal ölçüm standartlarına tabiidir. Bu kalibrasyon certificate documents the measurability to national standards, which realize the unit of measurement according to the International System of Units (SI). Kalibrasyon laboratuvarı olarak faaliyet gösteren Netes Mühendislik ve Tic. A.Ş. Kalibrasyon Laboratuvarı, Türk Standartları Enstitüsü (TSE) tarafından AB-0030-K akreditasyon dıģında numarası ile TS EN ISO/IEC 17025:2017 standardına göre akredite edilmiştir. Jener. Akreditörler ve Dış Tic. A.Ş. Calibration Laboratory accredited by TÜRKAK under registration number AB-0030-K for TS EN ISO/IEC 17025:2017 as Calibration Laboratory. Türk Akreditasyon Kurumu (TÜRKAK) kalibrasyon sertifikalarının tanımlı olduğu konularda Avrupa Akreditasyon Birliği (EA) ile Çok Taraflı Anlaşma ve Uyumlaşma Laboratuvarı Akreditasyon Birliği (ILAC) ile benzerlik tanıma anlaşmaları imzalamıştır. Turkish Accreditation Agency (TÜRKAK) is a signatory to the European co-operation for accreditation (EA) Multilateral Agreement (MLA) and to the International Laboratory accreditation Cooperation (ILAC) Mutual Recognition Arrangement (MRA) for the Recognition of calibration certificates. Ölçüm sonuçları, genişletilmiş ölçüm belirsizlikleri ve kalibrasyon miktarları bu sertifikasyon belgesinde ayrı ayrı belirtilen sayfalarca verilmektedir. The measurement, the uncertainty with confidence probability and calibration methods are given on the following pages which are part of this certificate. Bu sertifika, laboratuvarın yazılı izni olmadan kopyalanıp çoğaltılamaz. İmzasız sertifikalar geçerli değildir. This certificate shall not be reproduced other than in full except with the permission of the laboratory. Calibration certificates without signatures are not valid.

Mühür/Kayıt Seal	Yayınlandığı Tarih Date	Kalibrasyonu Yapan Calibrated by	Onaylayan / Tarih Approved / Date
	27.5.2020	Yılmaz ALÇIN	Ozan FİLİZ 27.5.2020

Bu sertifika laboratuvarın yazılı izni olmadan çoğaltılamaz. İmzasız ve mühürsüz sertifikalar geçerli değildir.
Tel: (0216) 340 50 50 pax Fax: (0216) 340 91 51 E-posta: www.neteskalibrasyon.com.tr E-posta: kalibrasyon@netes.com.tr
KİT FORSİSREV/01/0019 Sayfa 107 Page NUMBER 173

Declaration of Conformity



design mould & consulting company

Seyrinhisar Mah. Demirciler Sit. 9.Cd. No:26
Zeytinburnu / İstanbul / Turkey
Phone : +90 212 416 65 05 - 08
info@alpress.com.tr - www.alpress.com.tr

Declaration of Conformity

Manufacturer
Alpress Kalıpcılık Danışmanlık İth. Ihr. Ltd. Stl.
Demirciler Sit. 9.cd no26 Zeytinburnu İstanbul, Turkey

Herewith declare that:


Products description: ICU Ventilator
Product Name: IDEA Aero Duo
Model Number: DV0104

is in conformity with the following standards:

- IEC 60601-1: 2012: Medical Electrical Equipment – Part 1,
- IEC 60601-1-2: 2014: Medical Electrical Equipment Part 1-2,
- IEC 60601-1-11: 2015: Medical Electrical Equipment Part 1-11,
- IEC 62304: 2015: Medical Device Software – Software Life Cycle Processes
- ISO 10993: Fifth Edition 2018-08: Biological Evaluation of Medical Devices - Part 1
- ISO 18562-1 First Edition 2017-03: Biocompatibility Evaluation of Breathing Gas Pathways in Healthcare Applications - Part 1
- ISO 18562-2 First Edition 2017-03: Biocompatibility Evaluation of Breathing Gas Pathways in Healthcare Applications - Part 2
- ISO 18562-3 First Edition 2017: Biocompatibility Evaluation of Breathing Gas Pathways in Healthcare Applications - Part 3
- ISO 18562-4 First Edition 2017-03: Biocompatibility Evaluation of Breathing Gas Pathways in Healthcare Applications - Part 4

Issued by: Alpress Kalıpcılık Danışmanlık İth. Ihr. Ltd. Stl.
Date: 20/06/2020
Place: İstanbul / Turkey

Murat Degirmencioglu
Dept. Manager
Resource and Development Dept.
Alpress Company


KALIPÇILIK DANIŞMANLIK İTH. İHR. LTD. STL.
Demirciler Sit. 9.cd no26 Zeytinburnu İstanbul / Turkey
Tel: +90 212 416 65 05 - 08
Fax: +90 212 416 65 05 - 08
info@alpress.com.tr - www.alpress.com.tr

Declaration of Conformity



design mould & consulting company

Seyrinhisar Mah. Demirciler Sit. 9.Cd. No:26
Zeytinburnu / İstanbul / Turkey
Phone : +90 212 416 65 05 - 08
info@alpress.com.tr - www.alpress.com.tr

Declaration of Conformity

Manufacturer
Alpress Kalıpcılık Danışmanlık İth. Ihr. Ltd. Stl.
Demirciler Sit. 9.cd no26 Zeytinburnu İstanbul, Turkey

Herewith declare that:

Products description: ICU Ventilator
Product Name: IDEA Aero Duo
Model Number: DV0104

is in conformity with the following standard:

- ISO 80601-2-12 First Edition 2011-04-15: Medical Electrical Equipment - Part 2-12: Particular Requirements for the Safety of Lung Ventilators - Critical Care Ventilators

Issued by: Alpress Kalıpcılık Danışmanlık İth. Ihr. Ltd. Stl.
Date: 20/06/2020
Place: İstanbul / Turkey

Murat Degirmencioglu
Dept. Manager
Resource and Development Dept.
Alpress Company


KALIPÇILIK DANIŞMANLIK İTH. İHR. LTD. STL.
Demirciler Sit. 9.cd no26 Zeytinburnu İstanbul / Turkey
Tel: +90 212 416 65 05 - 08
Fax: +90 212 416 65 05 - 08
info@alpress.com.tr - www.alpress.com.tr



Thank You

Eon Group (Thailand) Co., Ltd. 

+66 2899 2400 

info@th.eon-group.com 

www.eon-group.com 

